

呼气末CO2传感器



Sidestream CO2 Module



Mainstream CO2 Sensor

Product Introduction

产品介绍

MED-LINKET's ETCO2 module is the ideal capnography solution. It is plug-in and play. The sensor features advanced technology for measuring instantaneous CO2 concentration, respiration rate, end-tidal CO2 concentration, inspired CO2 concentration. It resists the vapour better with patented technology, making the measured result more accurate.

美的连ETCO2探头是临床监测呼吸二氧化碳的完美解决方案，该产品即插即测，采用先进的非分光红外技术，可以测量被测对象的瞬时CO2浓度，呼吸率，呼气末CO2值，吸入CO2浓度。采用专利除水技术，更好的减小水汽干扰，使测量结果更精准。

Application Field

应用场合

- Monitoring the patient's respiratory status
监测患者的呼吸状况
- Assists in determining when to intubate or extubate
帮助医生判断何时通上呼吸机，何时撤离呼吸机
- Verification of ET tube placement
检验确认ET管位置是否放置正确
- Alerts if accidental extubation occurs
当发生意外拔管时可以及时报警
- Ventilator disconnect detection
帮助确认呼吸机通气环路意外脱落
- Verification of ventilation during transport
在患者转移过程中及时确认呼吸通路是否正常

Technical Specifications

性能指标

Transducer Type 传感器类型	Sidestream CO2 module(UN-E800) and Mainstream CO2 sensor(UN-E900) UN-E800 旁流呼气末CO2模块, UN-E900 主流呼气末CO2探头	
Principle of Operation 测量原理	Non-dispersive infrared (NDIR) 非分光红外技术, 测量反应速度快	
Initialization Time 初始化时间	Capnogram displayed in less than 20 seconds. 上电20秒内即会显示测量波形 Reaching full specifications within 2 minutes. 上电2分钟内达到最好精度	
Co2 Measure Range CO2 测量量程	0-150mmHg, 0-19.7%, 0-20kPa	
Co2 Resolution Measure CO2测量分辨率	0.1mmHg	
Co2 Accuracy CO2 测量精度	0 - 40 mmHg 41 - 70 mmHg 71 - 100 mmHg 101 - 150 mmHg	±2 mmHg ±5% of reading ±8% of reading ±10% of reading
Respiratory Rate 呼吸率	Range测量范围 Accuracy测量精度	3-150 BPM ±1 breath
Calibration 校准要求	No routine user calibration required, The mainstream sensor an airway adapter zero is required when changing to a different airway adapter 不需要用户频繁校准。对于主流CO2探头, 更换适配器后需要进行一次校零操作。	
Compensation 补偿参数	Barometer Pressure 450-850 mmHg. Operator selectable O2, N2O, He compensation. 压力补偿范围: 450-850mmHg, 可选择的O2, N2O, He气体补偿。	
Power Requirements 电源要求	5.0 VDC ±5%, power consumption less than 1.3W typical	
Temperature / Humidity 温度与湿度要求	Operating:0°C to 40°C, 10 to 90% RH 工作环境:温度要求0-40°C, 湿度要求10%-90% RH Storage:-40°C to 70°C, <90% RH 存储环境:温度要求-40°C-70°C, 湿度要求<90% RH	
Data Interface 数据接口	Rs232, 19200 baud rate, standard N-8-1. Fully compatible with Respironics CO2 module's communication protocol. 兼容伟康8PIN LEMO接口及其软件通信协议。	
Regulatory 标准	Designed to comply YY0601 and ISO21647. 所有测量值的精度符合YY0601和ISO21647的要求。	